




QUASAR™
MLC Beam
Geometry Phantom



QUASAR™
Respiratory
Motion QA



QUASAR™
Multi-Purpose
Body Phantom



QUASAR™
Verification Phantom



QUASAR™
Penta-Guide Phantom

QUASAR™ Penta-Guide Phantom

[Specifications](#) • [Brochure Download](#) • [References](#)

Specifications

- Cube, 16 cm, acrylic, 5 kg
- 4 x 4, 10 x 10, and 12 x 12 cm² light field alignment
- Laser alignment lines
- Built-in bubble level
- Excel spreadsheet for recording, analyzing and tracking results
- Internal imaging and registration markers designed to minimize CT artifacts
- User's Guide and cardboard container for storage and handling are included



Ordering Information

100-1009 QUASAR Penta-Guide Phantom